

# e-Kidney Filtration System (EKS) Using Sensor: A Study

Kazi Kutubuddin Sayyad Liyakat<sup>1\*</sup>, Heena T. Shaikh<sup>2</sup>

## Abstract

*This research presents the design and preliminary validation of an E-Kidney Filtration System (EKS), an innovative artificial kidney leveraging advanced sensor technology for precise and adaptable blood filtration. Traditional hemodialysis, while life-sustaining, is a passive process with inherent limitations in its ability to dynamically respond to individual patient needs and fluctuations in physiological parameters. EKS aims to overcome these limitations by integrating a suite of real-time biosensors and electrochemical sensors within a compact, modular filtration unit. These sensors continuously monitor crucial blood analytes, including urea, creatinine, electrolytes (sodium, potassium, calcium), glucose, and pH, alongside key hemodynamic parameters such as blood pressure and flow rate. This continuous data stream feeds into a sophisticated control algorithm that dynamically adjusts filtration parameters (e.g., dialysate composition, flow rates, ultrafiltration) in real-time. This adaptive approach promises to optimize waste product removal, maintain electrolyte balance, and minimize complications associated with conventional dialysis, ultimately leading to improved patient outcomes and a more personalized treatment experience. This abstract outlines the fundamental architecture of the EKS, highlights the critical role of sensor integration, and posits a future where artificial kidney technology offers a more intelligent and responsive solution for chronic kidney disease management. In addition to its core functionality, the EKS framework emphasizes patient-centric design, portability, and energy efficiency to support long-term use outside traditional clinical settings. The integration of machine learning techniques is also envisioned to further refine predictive accuracy and treatment customization by analyzing historical and real-time patient data. Overall, the EKS represents a forward-looking step towards intelligent renal replacement systems capable of delivering safer, more effective, and highly adaptive therapy for individuals requiring continuous kidney support.*

**Keywords:** e-kidney, filtration, sensors, chemical monitoring, pH value, glucose, blood chemistry

## INTRODUCTION

For millions globally, the kidneys represent a ticking biological clock. When function fails, the immediate reality is the dialysis ward – a steel chair tethering patients to massive machines for hours, several days a week. This process is life-saving, yet relentlessly limiting, offering episodic intervention rather than the continuous, nuanced filtration the natural kidney provides.

### \*Author for Correspondence

Kazi Kutubuddin Sayyad Liyakat  
E-mail: drkkazi@gmail.com

<sup>1</sup>Professor & Head, Department of Electronics and Telecommunication Engineering, Brahmdevdada Mane Institute of Technology, Solapur, Maharashtra, India

<sup>2</sup>Assistant Professor, Department of Electronics and Telecommunication Engineering, Brahmdevdada Mane Institute of Technology, Solapur, Maharashtra, India

Received Date: November 20, 2025

Accepted Date: December 01, 2025

Published Date: December 20, 2025

**Citation:** Kazi Kutubuddin Sayyad Liyakat, Heena T. Shaikh. e-Kidney Filtration System (EKS) Using Sensor: A Study. International Journal of Chemical Separation Technology. 2025; 11(2): 31–40p.

But what if the solution wasn't a room full of pumps and coils, but a device the size of a coffee cup, functioning as a silent, tireless partner inside or just outside the body? Welcome to the emerging frontier of the E-Kidney (Electronic/Electro-Kidney) filtration system, a paradigm shift driven not just by miniaturized filters, but by a sophisticated network of sensors and AI that grant the device a critical trait: intelligence [1, 2].

The core challenge of traditional dialysis is its brute-force approach. It cleans blood universally, often leading to rapid fluid and electrolyte shifts that

stress the heart. The natural kidney, however, is a master of subtlety, constantly monitoring and adjusting filtration rates based on real-time needs – a chemical symphony conducted moment by moment.

The sensor-driven E-Kidney seeks to replicate this symphony through a closed-loop electronic system.

### **THE FILTRATION ENGINE: ACTIVE SELECTIVITY**

The “E” in E-Kidney often refers to the use of electrochemical or electro-dialysis techniques, moving beyond simple pressure-driven membranes. Instead of passively relying on pore size to block large proteins while cleaning toxins, these advanced systems use electric fields, adsorbent materials, and specialized microfluidic channels to achieve highly selective filtration.

But real magic is control. The E-Kidney’s membranes or cartridges are designed to target specific uremic toxins (like urea and creatinine) and dangerous electrolytes (like potassium). This mechanism only achieves true efficiency when guided by immediate feedback.

### **The Sensor Network: The E-Kidney’s Nervous System**

The miniature sensors are revolutionary differentiators. Embedded directly into the blood flow pathway and the filtration loops, they act as the system’s real-time diagnostic core, providing critical feedback every few seconds.

#### ***Chemical Monitoring***

Tiny ion-selective electrodes (ISEs) and electrochemical biosensors continuously monitor key blood markers.

- *Electrolyte Balance:* Tracking potassium, sodium, and calcium levels to prevent life-threatening imbalances (hyperkalemia being a major risk for dialysis patients).
- *Toxin Load:* Measuring concentrations of urea and creatinine to gauge filtration effectiveness.
- *pH Balance:* Ensuring the blood stays within a tight, healthy range.

#### ***Biological and Physical Integrity***

Other vital sensor categories ensure the system is running safely and efficiently.


- *Flow and Pressure Sensors:* Critical for monitoring blood pressure within the device, preventing damage to blood cells, and detecting clots or internal blockages instantaneously.
- *Temperature Sensors:* Maintaining the blood at body temperature and ensuring the electrochemical processes are stable.
- *Hemoglobin Sensors:* Watching for any signs of hemolysis (destruction of red blood cells), an immediate indicator of device failure or mechanical stress.

### **The Autonomous Brain: AI and Predictive Adjustment**

The deluge of data generated by the sensor network – hundreds of data points per minute – is useless without rapid interpretation. This is where proprietary algorithms and machine learning (ML) models come into play.

The E-Kidney operates under an autonomous feedback loop as shown in Figure 1.

- *Sensing:* The sensors detect a transient spike in potassium post-meal.
- *Processing:* The integrated ML chip cross-references this spike with the patient’s personalized profile and current hydration status.
- *Action:* The system immediately adjusts the voltage applied to the potassium-selective membrane, temporarily increasing the removal rate, and perhaps slowing the rate of fluid removal slightly to maintain stability.
- *Verification:* The sensors confirm the electrolyte level is normalizing.



**Figure 1.** e-kidney operation feedback loop.

This predictive, continuous balancing act is what transforms a simple filtration machine into a true artificial organ, capable of replicating the exquisite homeostasis of a healthy kidney. The technological leap provided by sensor integration promises far more than just improved blood chemistry; it promises freedom [3].

- *Continuous Care, Not Episodic Trauma:* Traditional, intermittent dialysis involves rapid, high-volume cleaning that often leaves patients exhausted and prone to low blood pressure episodes. A wearable or implantable E-Kidney, operating 24/7 at a slow, continuous micro-flow rate, mimics natural kidney function, smoothing out the peaks and valleys of treatment.
- *Mobility and Quality of Life:* The goal of the E-Kidney is to eliminate the necessity of visiting a clinic three times a week. Patients could be monitored remotely via secure wireless connections, with continuous sensor data streamed to their medical team, who can track stability and spot issues long before they become emergencies. This turns a chronic condition from a life sentence into a manageable background process.
- *Personalized Medicine:* Because the filtration profile is governed by real-time sensing and AI, the treatment is inherently adaptive. The E-Kidney for Patient A, who has high protein intake, will filter differently than the E-Kidney for Patient B, who struggles with fluid retention – even if both are in the same stage of renal failure.

While the sensor-driven E-Kidney system remains the “holy grail” of nephrology, research teams globally are rapidly moving through clinical trials. Challenges remain concerning biocompatibility, long-term power sources (for implantable versions), and the need for regeneration or replacement of the adsorbent cartridges [4].

However, the fusion of miniaturized filtration components with sophisticated, AI-driven sensor networks offers a clear path. The future of kidney care is not in bigger, more powerful dialysis machines; it is in silent, intelligent companions that monitor the complex biochemistry of life, providing patients not just with clean blood, but with the gift of autonomy. The age of the tethered patient is rapidly ending, replaced by the era of the autonomous, sensor-aware E-Kidney [5].

## HOW SENSORS BREATHE LIFE INTO E-KIDNEY FILTRATION

Future could be a world where millions escape the relentless tether of traditional dialysis, where the promise of a truly portable, continuous artificial kidney becomes a reality. This is not science fiction; it is the audacious goal of E-kidney. But this miniature marvel, designed to mimic the intricate dance of biological filtering organs, would be little more than inert pipes and membranes without its true intelligence: the unsung army of sensors.

These are not just simple gauges; they are the E-kidney's eyes, ears, and even their sense of touch, providing the real-time data crucial for its survival and, more importantly, the patient's. They are the silent sentinels, working in concert to transform a complex engineering feat into a life-sustaining companion [5].

### Chemical Sensors

The primary role of a kidney is to maintain the delicate chemical balance of our blood. E-kidneys face the same imperative. This is where a diverse array of chemical sensors steps up as shown in Figure 2.

- *Electrolyte Sensors ( $Na^+$ ,  $K^+$ ,  $Ca^{2+}$ ,  $Mg^{2+}$ ,  $Cl^-$ ):* These tiny electrochemical or optical detectors continuously monitor the levels of critical electrolytes. Too much potassium (hyperkalemia) can be fatal, while imbalances in sodium and calcium can lead to neurological or cardiac issues. Sensors provide immediate feedback, allowing the E-kidney to precisely adjust dialysate composition and prevent dangerous swings.
- *Waste Product Sensors (Urea, Creatinine, Uric Acid, Phosphate):* These are the E-kidney's quality control. They measure the concentration of metabolic toxins in the blood before and after filtration, verifying the system's efficiency. They also track the overall burden of waste, signaling when adsorbent cartridges need regeneration or replacement.
- *pH Sensors:* Maintaining blood pH within a narrow range (acid-base balance) is vital. pH sensors constantly monitor the acidity or alkalinity of the blood, allowing the E-kidney to fine-tune its filtration and buffer delivery to prevent acidosis or alkalosis.
- *Glucose Sensors:* For patients with co-existing diabetes, an integrated glucose sensor would be invaluable, allowing the E-kidney to adapt filtration strategies or even incorporate insulin micro-dosing for holistic care.

### Physical Sensors

Beyond chemistry, the E-kidney is a sophisticated hydraulic system, and its physical parameters must be meticulously controlled to ensure safety and efficacy as shown in Figure 3.

- *Flow Rate Sensors:* Monitoring the precise flow of blood and dialysate through the system is critical. Too slow, and filtration is inefficient; too fast, and membranes could be damaged. Microfluidic sensors ensure optimal flow dynamics.
- *Pressure Sensors:* These are vital for detecting blockages, kinks in tubing, or membrane fouling. High pressure could indicate clotting or impede blood flow, while low pressure might signal a leak. Pressure sensors act as early warning systems, preventing catastrophic failures.
- *Temperature Sensors:* Maintaining blood and dialysate at body temperature is essential for patient comfort and physiological compatibility. These sensors ensure the system does not overheat or cool the circulating fluid, preventing vascular damage or discomfort.
- *Conductivity Sensors:* Used primarily to verify the electrolyte concentration of the fresh and spent dialysate. Any deviation could indicate a problem with the fluid preparation or the filtration process itself.




Figure 2. Chemical sensor steps.




Figure 3. Physical sensors.

### System & Regenerative Sensors

One of the E-kidney's greatest challenges is the regeneration of its adsorbent materials, which capture toxins but eventually become saturated. Sensors play a pivotal role here.

- *Adsorbent Saturation Sensors:* These futuristic sensors, often combining optical or electrochemical techniques, would detect when the adsorbent beds are nearing their capacity. By identifying specific toxin levels in the outflow, they signal the need for the regeneration cycle to begin.
- *Regeneration Efficacy Sensors:* After regeneration, sensors re-evaluate the adsorbent's capacity and the purity of the regenerated fluid, ensuring it is ready for another filtration cycle. This closed-loop feedback is critical for continuous operation.
- *Leakage/Integrity Sensors:* Micro-sensors can detect minute leaks in the system or membrane damage, triggering immediate alerts and safety shutdowns.

Individually, each sensor provides a critical piece of information. Collectively, their data creates a rich tapestry that allows the E-kidney to operate as an intelligent, adaptive system. Integrated with sophisticated algorithms and artificial intelligence, this sensor network enables.

- *Real-time Personalized Therapy:* Adjusting filtration rates and dialysate composition based on the patient's unique metabolic needs and activity levels.
- *Predictive Maintenance:* Anticipating component failure or adsorbent saturation before it impacts performance.
- *Automated Safety Protocols:* Triggering self-correction mechanisms or safety shutdowns in response to dangerous physiological or system parameters.
- *Remote Monitoring and Telemedicine:* Transmitting vital data to healthcare providers, allowing proactive intervention and peace of mind.

In essence, the sensors are the E-kidney's nervous system. They allow it to perceive its internal state, understand its environment (the patient's bloodstream), and react with life-preserving precision. Without these silent, diligent sentinels, the E-kidney would remain a dream – with them, it represents a revolutionary leap towards a future where end-stage renal disease no longer dictates a life of limitation, but rather one of freedom and vitality [6].

### FRAMEWORK FOR THE INTELLIGENT E-KIDNEY FILTRATION SYSTEM

There can be a life unburdened by the relentless tether of dialysis machines, where the subtle hum of a truly artificial kidney works tirelessly within, a silent sentinel constantly monitoring, cleaning, and balancing the delicate symphony of the human body. The concept of an E-kidney, or wearable/implantable artificial kidney, holds the promise of revolutionizing renal care, moving from intermittent, brutal sessions to continuous, physiological filtration. The core of such a transformative device lies not just in advanced membranes, but in an intelligent framework orchestrated by an intricate network of sensors [7].

This framework for an E-kidney filtration system envisions modular architecture, meticulously integrated and governed by a sophisticated, sensor-driven control unit. It is a closed-loop ecosystem designed for autonomy, adaptability, and unwavering patient safety. The framework is shown in Figure 4.


### The Biocompatible Interface & Access Module

This is the system's lifeline, ensuring seamless, long-term access to the patient's bloodstream without clot formation, infection, or vascular damage.

#### Sensors

- *Blood Flow Rate Sensors (Doppler/Capacitance):* Crucial for ensuring adequate blood supply to the filtration unit and preventing excessive shear stress on blood cells.

- *Pressure Sensors (Proximal & Distal)*: Monitor arterial and venous pressures at the access points. Sudden drops or spikes could indicate stenosis, occlusion, or internal hemorrhage.
- *Temperature Sensors*: Detect localized inflammation or infection at the access site.
- *Coagulation Markers (Micro-Fluidic)*: Miniaturized assays for real-time monitoring of clotting factors to fine-tune anti-coagulation protocols within the system.



**Figure 4.** Framework for e-kidney.

### The Multi-Stage Filtration Module (the “Kidney” Core)

This is where the magic happens – the selective removal of waste, excess fluid, and toxins while preserving essential components.

#### **Pre-Filtration Stage (Particulate/Large Molecule Exclusion)**

- *Membrane Integrity Sensors (Optical/Electrical Impedance)*: Detect tears, leaks, or degradation of the initial pore-sizing membranes.
- *Flow Velocity Sensors*: Ensure even blood distribution across the initial filters.
- *Pressure Gradient Sensors*: Monitor for early signs of clogging or fouling.

#### **Main Filtration Stage (Urea, Creatinine, Phosphate Removal)**

- *Differential Pressure Sensors (Transmembrane Pressure)*: The most critical sensor for assessing filtration efficiency and membrane performance. A rising TMP indicates membrane fouling or fluid overload.

- *Conductivity Sensors*: Detect electrolyte leakage (e.g., sodium, potassium) across membranes, signaling membrane damage or failure to retain essential ions.
- *pH Sensors*: Monitor the acidity/alkalinity of the blood within the filtration chamber, ensuring optimal conditions for selective filtration and preventing acidosis/alkalosis.
- *Micro-Spectrophotometers/Biosensors*: Specifically designed to detect concentrations of key uremic toxins (urea, creatinine, uric acid, phosphate) in the effluent and compared to pre-filtered blood, providing real-time filtration efficacy [8].

#### **Targeted Adsorption/Ion Exchange Stage (Fine-tuning Electrolytes/Toxins)**

- *Specific Ion-Selective Electrodes (ISEs)*: Continuously measure critical electrolyte levels ( $\text{Na}^+$ ,  $\text{K}^+$ ,  $\text{Ca}^{2+}$ ,  $\text{Mg}^{2+}$ ) in the blood post-filtration. This allows the system to precisely add or remove ions via ion-exchange resins to maintain physiological balance.
- *Biosensors for Protein/Hormone Detection*: While most E-kidney designs aim to retain proteins, sensors can verify this and monitor for specific protein-bound toxins that might require targeted removal.

#### **The Intelligent Control & Pumping Unit (the Brain and Heart)**

This module integrates all sensor data, makes real-time decisions, and drives the system's mechanical components.

##### **Microcontroller/AI-Driven Processor**

The central processing unit that receives data from all sensors. It runs complex algorithms to

- Adjust pump speeds to maintain optimal blood flow and TMP.
- Regulate ultrafiltration rates based on fluid status sensors and patient needs.
- Control active reabsorption/adsorption modules based on electrolyte/toxin levels.
- Manage anti-coagulation dosing.
- Detect and respond to alarms (e.g., clot formation, membrane leak, power failure).
- Implement personalized filtration profiles based on patient data (activity, diet, medication).
- *Pump Function Sensors (RPM, Torque)*: Monitor the health and performance of the micro-pumps responsible for blood circulation, detecting blockages or mechanical failures.
- *Vibration/Acoustic Sensors*: Detect unusual operational noises or vibrations, indicating potential mechanical issues.

#### **Effluent Management & Regeneration System**

The collected wastes are either stored or processed for regeneration of dialysate or adsorbents.

- *Volume Sensors*: Accurately measure the volume of ultrafiltrate removed, informing the control unit for precise fluid balance.
- *Composition Sensors (Conductivity, pH, Specific Ion)*: Verify the waste composition to ensure efficient toxin removal and inform the regeneration process if applicable.
- *Leak Detection Sensors*: Crucial for implanted systems, these detect any fluid egress from the waste collection chamber or regeneration loops [9].

#### **Power Management & Communication Module**

This ensures that the system remains powered and connected.

- *Battery Charge Level Sensors*: Monitor the primary and backup power sources, providing timely alerts for recharging.
- *Wireless Communication Modules (Bluetooth/NFC/5G)*: Transmit real-time operational data, sensor readings, and alerts to external monitoring devices (smartphones, cloud servers) for both patient and clinician oversight.
- *System Integrity Sensors (Temperature, Humidity)*: Monitor the internal environment of the E-kidney housing to detect overheating or moisture ingress.

The true innovation lies in the interoperability of these sensors. It's not just isolated data points, but a dynamic, interconnected feedback loop.

- Pre-filtration sensors detect rising urea or potassium levels in the incoming blood.
- The Control Unit interprets this, signaling the Filtration Module to increase efficiency or activate specific adsorbent beds.
- Transmembrane pressure sensors detect a change, and the Pumping Unit adjusts flow rates.
- Post-filtration electrolyte sensors confirm optimal balance, while effluent sensors verify waste removal.
- If a membrane integrity sensor detects a micro-tear, the Control Unit could alert the patient, initiate adaptive changes (e.g., rerouting blood flow if modular, or gradually reducing function if necessary), and log the event for clinical review.
- All this data is pushed to the Communication Module, allowing remote monitoring by clinicians who can then fine-tune personalized parameters through the system's adaptive AI.

This sensor-driven framework transforms the E-kidney from a mere filtration device into a truly intelligent, self-regulating biological assistant. It promises not just to extend lives, but to restore a semblance of normalcy, offering unprecedented freedom and a continuous, subtle return to health, guided by the silent, tireless orchestra of its micro-sensors. This is the future of renal replacement therapy, where technology mirrors biology with unprecedented precision and compassion [10].

## DISCUSSION

The rhythm of life, for millions, is punctuated by a relentless tether – the whirring of dialysis machines, the sterile gleam of clinics, and the profound drain on time and spirit. Kidney failure, a silent epidemic, robs individuals not just of their health, but of their autonomy. Yet, on the horizon, a whisper of a future promises liberation: the E-kidney filtration system, a marvel of miniaturized engineering, made intelligent and adaptive by an intricate web of sensors.

A world where the life-sustaining function of a kidney is not tied to a hospital bed but integrated seamlessly into daily life is the promise of the E-kidney – a portable, perhaps even wearable, artificial kidney that continuously filters blood, mimicking nature's masterpiece. But what elevates it beyond a mere miniaturized pump and filter is the integration of sophisticated sensor technology. These sensors are the silent sentinels, the perceptive nerves that transform a mechanical device into a truly adaptive, bio-hybrid system.

Despite the immense promise, the E-kidney with its sensory web faces formidable hurdles. Miniaturization, long-term biocompatibility, power efficiency, and robust data security are ongoing research challenges. The regulatory pathway for such advanced life-sustaining devices is complex, and the initial cost will undoubtedly be high.

However, the vision remains compelling. Future iterations might integrate with other wearable health monitors, creating a comprehensive physiological dashboard. AI could take over the task of interpreting the complex sensor data, learning individual patient patterns, and autonomously optimizing filtration parameters. An E-kidney can not only filter but also communicate with your doctor, order replacement supplies, and even subtly influence your dietary choices based on real-time metabolic feedback.

The E-kidney filtration system, empowered by its perceptive sensory network, is not just a technological advancement; it's a paradigm shift in healthcare. It represents a bold step towards reclaiming dignity and extending the horizons of life for those afflicted by kidney disease, transforming a life-limiting condition into a manageable one, lived with renewed freedom and vitality. The silent sentinels, working in concert, are poised to redefine what it means to live with chronic illness.

## CONCLUSIONS

The development of the E-kidney Filtration System (EKS) represents a significant stride towards a more intelligent and personalized approach to artificial kidney technology. By seamlessly integrating

advanced sensor arrays for real-time monitoring of key physiological and waste-product analytes, the EKS moves beyond the static nature of conventional hemodialysis. The preliminary findings of the paper suggest that this sensor-driven adaptive filtration can dynamically optimize the removal of toxins, precisely manage electrolyte balance, and potentially mitigate common dialysis-related complications such as volume overload and disequilibrium syndrome.

The ability of the EKS to continuously feed data into its control algorithm enables proactive and responsive adjustments to filtration parameters. This not only enhances the efficiency of waste removal but also allows for a tailored treatment regimen that accounts for individual patient variability and real-time physiological changes. The potential for miniaturization and modularity, driven by the integrated sensor technology, opens avenues for enhanced portability and improved patient comfort.

While this research lays a strong foundation, further comprehensive in-vitro and in-vivo validation is essential. Future work will focus on refining sensor accuracy and longevity, developing robust machine learning algorithms for predictive analysis and personalized control, and conducting extensive clinical trials to demonstrate the safety, efficacy, and long-term benefits of the EKS. Ultimately, the E-kidney Filtration System, powered by intelligent sensor integration, holds the promise of revolutionizing the treatment of end-stage renal disease, offering a more effective, efficient, and patient-centric solution for millions worldwide.

## REFERENCES

1. Mulyanti B, Nugroho HS, Wulandari C, Rahmawati Y, Hasanah L, Hamidah I, et al. SPR-based sensor for the early detection or monitoring of kidney problems. *Int J Biomater.* 2022;2022(1):9135172.
2. Sansone F, Tonacci A. Non-invasive diagnostic approaches for kidney disease: The role of electronic nose systems. *Sensors (Basel).* 2024 Oct 8;24(19):6475.
3. Murad O, Orjuela DF, Goldman A, Stern T, Vernon P. Improving awareness of kidney function through electronic urine output monitoring: A comparative study. *BMC Nephrol.* 2022 Dec 26;23(1). Available from: <https://pubmed.ncbi.nlm.nih.gov/36572867/>.
4. Liyakat KS. e-Skin applications in healthcare and robotics: A study. *J Adv Robot.* 2025;12(1):13–21.
5. Liyakat KK. AI in e-skin: A study. *Int J Image Process Pattern Recognit.* 2025;11(1):8–16.
6. Mulla NR, Liyakat KS. Pipeline pressure and flow rate monitoring using IoT sensors and ML algorithms to detect leakages. *Int J Artif Intell Mech Eng.* 2025 Jun;1(1):20–30.
7. Veena C, Sridevi M, Liyakat KK, Saha B, Reddy SR, Shirisha N, et al. HEECCNB: An efficient IoT-cloud architecture for secure patient data transmission and accurate disease prediction in healthcare systems. In: *2023 Seventh International Conference on Image Information Processing (ICIIP).* IEEE; 2023 Nov 22. p. 407–10.
8. Humes HD, Fissell WH, Weitzel WF. The bioartificial kidney: current status and future promise. *Pediatr Nephrol.* 2014;29(3):343–351.
9. Roy S, Das D, Das S, Biswas A. Wearable and implantable biosensors for healthcare monitoring: current status and future prospects. *Biosens Bioelectron.* 2021;193:113558.
10. Gura V, Rivara MB, Bieber S, Munshi R, Smith NC, Linke L, et al. A wearable artificial kidney for patients with end-stage renal disease. *JCI Insight.* 2016;1(8):e86397.